

Evaluation of VIDAS® *C.difficile* Toxin A IgG and VIDAS® *C.difficile* Toxin B IgG prototypes on blood bank sera characterized with cell cytotoxicity neutralisation assay

D. Rolland, M.C. Bastide, M.L. Coutin, C. Compagnion, F. Cottin, M. Dupin, A. Foussadier
bioMérieux, Marcy l'Etoile, France

INTRODUCTION

Clostridium difficile (CD) is a causative agent of antibiotic-associated pseudo membranous colitis (PMC). Since 2000, there has been a striking increase in the prevalence of *Clostridium difficile* infection and in associated mortality in the United States, Canada, and Europe. In this context, bioMérieux commercialize currently a VIDAS® *C.difficile* Toxin A&B test for toxins A and B detection. It seems necessary to investigate the impact of antibody responses to *Clostridium difficile* toxins on the risks of colonization, diarrhea and asymptomatic carriage.

OBJECTIVES

The aim of this study was to evaluate the results obtained on blood bank samples characterized with cell cytotoxicity neutralisation assay (CCNA) on the VIDAS® *C.difficile* Toxin B IgG (CDB IgG) prototype and on the VIDAS® *C.difficile* Toxin A IgG (CDA IgG) prototype for the determination of *C. difficile* serological status.


MATERIALS AND METHODS

VIDAS test:

The prototypes VIDAS CDB IgG and VIDAS CDA IgG use specific native toxin B and toxin A antigens respectively coated on the solid phase with the help of specific monoclonal antibodies against these toxins respectively. During the first step of the test, the specific IgG are captured by this solid phase. The complexes are revealed during the second step using an anti-human IgG conjugated to Alkaline Phosphatase (PAL). The result is presented in signal (RFV = relative fluorescent value). The serum volume is 100 µl for each test. The results are obtained in 35 minutes. Results are given as a relative fluorescence value for each sample and compared to results obtained with cell cytotoxicity neutralisation results.

For this study, the cut off is established at 300 RFV.

VIDAS test format



Cell Cytotoxicity Neutralisation Assay (CCNA):

100 µl of human sera, diluted to 1/20 are dispensed in duplicate on Vero cell microtiter plates and then, ten fold logarithmic dilution of *C. difficile* toxin B (110 µg/ml) are dispensed, 100 µl per well. The first well of each row received the dilution medium as negative control. Dilution of toxin B are dispensed in two wells as positive control. Microtiter plates are incubated at 37°C under 5% CO₂ during 48 h. The cell monolayer are checked by microscopy after 24 and 48 h in order to note the typical cell rounding effect. The results are expressed as the lower dilution of toxin B neutralized by the specific antibodies present in the sera.

Samples:

Forty two sera from blood bank origin are tested on VIDAS with each prototype. These sera were characterized with the cell cytotoxicity neutralisation assay currently used in bioMérieux.

RESULTS

Of the 18 samples characterized as positive in the cell cytotoxicity neutralisation assay, 18 samples showed positive results with the VIDAS CDB IgG prototype and 9 with the VIDAS CDA IgG prototype.

Of 24 samples characterized as negative in the cell cytotoxicity neutralisation assay, 18 samples showed negative results with the VIDAS CDB IgG prototype and 17 with the VIDAS CDA IgG prototype.

The agreement with the cell cytotoxicity neutralisation assay is 86% for VIDAS CDB IgG prototype and 62% for VIDAS CDA IgG prototype on the 42 samples.

The results are detailed in the tables below.

Samples	CTA neutralisation 1/20	anti-toxin A IgG	anti-toxin B IgG
	cut off	signal positive > 300 RFV	
CD029		1303	825
CD032		280	201
CD033		570	124
CD038		903	137
CD031		298	40
CD039		1090	389
CD024		390	660
CD027		137	41
CD025		377	1088
CD036		424	861
CD038		123	503
CD023		161	40
CD028		30	11
CD028		188	251
CD027		50	58
CD038		17	11
CD029		150	40
CD025		238	58
CD028		47	7
CD029		65	77
CD021		330	385
CD040		194	113
CD04		440	217
CD032		108	24

Samples	CTA neutralisation 1/20	anti-toxin A IgG	anti-toxin B IgG
	cut off	signal positive > 300 RFV	
CD034	- 4 logs	1186	4944
CD036		1290	5239
CD043	- 3 logs	2259	2616
CD041		254	2251
CD038		1545	1967
CD044		136	5472
CD035		3494	4395
CD030		195	3946
CD045	- 2 logs	619	3157
CD034		581	2284
CD05		119	2248
CD042		159	1800
CD031		1247	609
CD030		34	1895
CD09	- 1 log	191	848
CD07		155	548
CD034		145	558
CD037		552	388

	n	Correlation	CI 95%
anti-toxin B IgG	36/42	85.71%	[71.46-94.57]
anti-toxin A IgG	26/42	61.9%	[45.64-76.43]

Sensitivity	anti-toxin B IgG	anti-toxin A IgG
Positive samples	18	9
Negative samples	0	9
%	100% [81.47-100]	50% [26.02-73.98]

Specificity	anti-toxin B IgG	anti-toxin A IgG
Negative samples	18	17
Positive samples	6	7
%	75% [53.29-90.23]	71% [48.91-87.38]

CONCLUSION

With the VIDAS *C. difficile* Toxin A IgG and the VIDAS *C. difficile* Toxin B IgG prototypes, it is possible to detect the *C. difficile* toxins A and B IgG in blood samples. The results obtained with the VIDAS *C. difficile* Toxin B IgG prototype are better correlated with the cytotoxicity neutralisation assay than the results with the VIDAS *C.difficile* Toxin A IgG prototype ($p<0.05$). The VIDAS *C. difficile* Toxin B IgG seems to be more relevant for studying the immunological response for *C. difficile* infection. Larger studies with different population in the context of *C. difficile* infection in hospital settings (patients at risk, colonized patients, infected patients and patients with asymptomatic carriage,...) will need to be conducted to establish the optimal neutralisation cut off of the VIDAS *C. difficile* Toxin B IgG.